


# CLEANING A PATIENT CHAIR


It is important to check the integrity of the cover, holes or rips can let in contaminated body fluids. Unzip the cover and check as per local policy and between patients as a minimum. If you find damage, you should report it.

1/6



Perform hand hygiene. Wear appropriate PPE.

2/6




Remove one wipe from Universal pack. Wipe from top to bottom, clean to dirty in a 'S' motion.

3/6

1. TOP SURFACE


2. FRAME

3. UNDERSIDE




Remove cushions, clean in order then replace any cleaned cushions.

4/6




Change wipe if it is dry or soiled and discard. Allow surfaces cleaned with Universal wipes to air dry.

5/6



Remove all PPE and discard as per local policy. Perform hand hygiene.

6/6



As per local policy, add cleaned date and sign the indicator tape/note.

© 2020 GAMA Healthcare. All rights reserved. No extract may be reproduced, stored or transmitted in any form or by any means without prior written permission.

Please ensure that you follow the equipment manufacturer's guidelines.

GHA200368